

*** Crude DNA sample preparation from tail biopsy**
Youxi Ai (July 19, 2006)

1. Put tails (~2mm) in 1.5ml tubes;
2. Add 300µl of 50mM NaOH and put in steamer (pre-heated 10') for 1 hr or 95°C heat block for 50 min;
3. Neutralize with 30µl of 1M Tris (pH 8.0) after cooled down and invert the tubes (**Do Not Vortex!**);
4. Quick spin for 5' at 13,000 rpm (then take 1.5µl sup. for PCR).

Reaction :

1.5 uL DNA
1.5 uL 10x PCR Buffer (NEB with MgCl₂)
3.0 uL 1.25mM dNTPs,
0.25 uL For Primer (10 uM)
0.25 uL Rev Primer (10 uM)
0.15 uL NEB Taq (part # M0273S, 5U/ul)
8.35 uL PCR water

total : 15 ul

PCR program,

1. 95C 5 min
2. 94C 1 min
3. 55C 30 sec
4. 72C 1 min
5. go back to step 2, repeat for 36 cycles
6. 72C 7 min
7. 12C forever

Primers

***HuR** (Loxp 385; WT 290; Δ 507):

F: 5'-AGGCAGATGAGCACATGTGA-3' (3455-84)

R: 5'-AGGCTCTGGGATGAAACCTA-3' (3821-40)

ΔR2: 5'-TACTGAGATGTTCTGGGAGG-3' (17657-17676)

SK1: (WT 300; KO 350) **annealing temp 60C**

F: 5'-AGAAGGCACTGGCTCCTCCAG-3'

R: 5'-TGTCACCCATGAACCTGCTGT-3'

Neo 5'-TCGTGCTTTACGGTATCGCCGCTCCC-3'

SK2 (KO 300; WT 700 ?)

F: 5'-GCACCCAGTGTGAATCGAGC-3'

R: 5'-TCTGGAGACGGGCTGCTTTA-3'

Neo 5'-CGCTATCAGGACATAGCGTT-3'

SPK1 f/f : (KO 600, flox 250, WT 120) **annealing temp 60C**

lox F: 5'-GGA CCT GGC TAT GGA ACC -3'

lox R1: 5' -ATG TTT CTT TCG AGT GAC CC -3'

lox R2: 5'- AAT GCC TAC TGC TTA CAA TAC C -3'

Edg-1: (WT ~400; lacZ 250)

F: 5'-CCATCCTCTGCAGGATCT-3'

R: 5'-TGCTGCGGCTAAATTCCATG-3'

Neo 5'-TCGCCTTCTTGACGAGTTCTTCTGAG-3'

S1P1/Loxp (Loxp 250, WT 200, deleted Δ/Δ 200):

loxP1: 5'- GAG CGG AGG AAG TTA AAA GTG - 3'

loxP2: 5' - CCT CCT AAG AGA TTG CAG CAA -3'

loxP3: 5' - GAT CCT AAG GCA ATG TCC TAG AAT GGG ACA -3'

S1P2^{-/-} (- 200; + 150):

For 5'-GCAGTGACAAAAGCTGCCGAATGCTGATG-3'

Rev 5'-AGATGGTGACCACGCAGAGCACGTAGTG-3'

Neo 5'-TGACCGCTTCCTCGTGCTTTACGGTATCG-3'

S1P3 : (KO 380; WT 130)

F: 5'-TCAGTATCTTCACCGCCATT-3'

R: 5'-AATTCACACTACGGTCCGCAGAA-3'

Neo 5'-GTGCAATCCATCTTGTTCAAT-3'

S1P2 f/f (fl/fl 444 ; WT 349; deleted Δ/Δ 216)

loxF1 5' CAA GCT CTC TAC CAA GTG AGT TAC 3' (3478-3501)

loxR1 5' TGACCCTTCATGTGACAGGCAGAA 3' (3921-3898)

loxR2 5' GGC CAA TTA CAA CCA CCA TAG A 3' (7284-7263)

ApoM : (WT 154 ; KO 154)

mApoM 51:CACCCAGCAACTCATCCTTT
mApoM 31: GCAGCCATGTTGAAGACAAT

Neo-51:GTAGCCGGATCAAGCGTATG
Neo-31:CTGTGCTCGACGTTGTCACT

ApoM Tg: (403)

hApoM-51: GGGACTTGAATTCCTCCACA
hApoM-31: TGAAGGGAGCACAGATCTCA

S1P1 f/s/f : (700)

209bF: 5'-CACTGCATTCTAGTTGTGGT
208bR: 5'- GTGTAGTTTCATCTTCAGCA

ApoE (WT 155 ; KO 245)

180 : 5' GCC TAG CCG AGG GAG AGC CG-3' Common
181 : 5' TGT GAC TTG GGA GCT CTG CAG C 3' WT Reverse
182 : 5' GCC GCC CCG ACT GCA TCT 3' Mutant Reverse

S5A (WT 400; KO 300)

S5a Fwd ; 5'-AAG TGC AGT GAG GCC AAA CA-3'
S5a Rev: 5'-GCA AGC ATG AGG CAA GRC TA-3'

gCre (300)

F: GATATCTCACGTA CTGACGG
R: TGACCAGAGTCATCCTTAGC


Rosa Cre (WT 550; Cre 350):

F: 5'-AAAGTCGCTCTGAGTTGTTAT-3'
R1: 5'-GCCAAGAGTTTGTCTCAACC-3'
R2: 5'-GGAGCGGGAGAAATGGATATG-3'

****LyzsM Cre** (WT 350; Cre 700) **annealing temp 56C**

For 5'-CCCAGAAATGCCAGATTACG-3'
Rev1 5'-CTTGGGCTGCCAGAATTTCTC-3
Rev2 5'-TTACAGTCGGCCAGGCTGAC-3'

ApoM KO and Tg genotyping



mApoM : WT

mApoM (in1-2)-51: CACCCAGCAACTCATCCTTT

mApoM (ex2 target)-31: GCAGCCATGTTGAAGACAAT

ApoM Neo

Neo-51: GTAGCCGGATCAAGCGTATG


Neo-31: CTGTGCTCGACGTTGCACT

Hu ApoM

hApoM-51: GGGACTTGAATTCCTCCACA

hApoM-31: TGAAGGGAGCACAGATCTCA

S1P receptors genotyping



S1P1 ^{f/f}

S1P1^{flox}-P1: 5' GAG CGG AGG AAG TTA AAA GTG 3'

S1P1 ^{flox} -P2: 5' CCT CCT AAG AGA TTG CAG CAA 3'

S1P1 ^{flox} -P3: 5' GAT CCT AAG GCA ATG TCC TAG AAT GGG ACA 3' (for deleted allele)

S1P2 ^{-/-}

For 5'-GCAGTGACAAAAGCTGCCGAATGCTGATG-3'

Rev 5'-AGATGGTGACCACGCAGAGCACGTAGTG-3'

Neo 5'-TGACCGCTTCTCTGCTTTACGGTATCG-3'

S1P3^{-/-}

S1P3 For : 5'-TCAGTATCTTCACCGCCATT-3'

S1P3 Rev: 5'-AATTCACTACGGTCCGCAGAA-3'


S1P3 Neo: 5'-GTGCAATCCATCTTGTTC AAT-3'

Generic Cre-ERT2

For: 5'-GAT ATC TCA CGT ACT GAC GG- 3'

Rev: 5'-TGA CCA GAG TCA TCC TTA GC-3'

ApoE, S5A and LysM Cre




ApoE: WT: 155 bp, KO: 245 bp
 S5A: WT: 400 bp, KO: 300 bp
 LysMCre: Mutant: 700 bp, WT: 300 bp

S5A
 S5a Fwd ; 5'-AAG TGC AGT GAG GCC AAA CA-3'
 S5a Rev: 5'-GCA AGC ATG AGG CAA GRC TA-3'

LysMCre
 For 5'-CCCAGAAATGCCAGATTACG-3'
 Rev1 5'-CTTGGGCTGCCAGAATTCTC-3'
 Rev2 5'-TTACAGTCGGCCAGGCTGAC-3'

ApoE
 180 : 5' GCC TAG CCG AGG GAG AGC CG-3' Common
 181 : 5' TGT GAC TTG GGA GCT CTG CAG C 3' WT Reverse
 182 : 5' GCC GCC CCG ACT GCA TCT 3' Mutant Reverse


S1P2^{f/f} genotyping



WT: 349 bp
f/f : 444 bp

S1P2_F1 5' CAA GCT CTC TAC CAA GTG AGT TAC 3' (3478-3501)
S1P2_R1 5' TGACCCTTCATGTGACAGGCAGAA 3' (3921-3898)
 Δ R2 5' GGC CAA TTA CAA CCA CCA TAG A 3' (7284-7263) $\Delta/\Delta=213$ bp


HuR genotyping



f/f : 385 bp
WT: 290 bp

F: 5'-AGGCAGATGAGCACATGTGA-3' (3455-84)
R: 5'-AGGCTCTGGGATGAAACCTA-3' (3821-40)
ΔR2: 5'-TACTGAGATGTTCTGGGAGG-3' (17657-17676)

S1P1^{f/s/f} genotyping



Target: 700bp

209bF: 5'-CACTGCATTCTAGTTGTGGT
208bR: 5'-GTGTAGTTTCATCTTCAGCA