

Motion Management in Kidney Diffusion-Weighted MRI

Simon K. Warfield, Ph.D.

Thorne Griscom Professor of Radiology

Department of Radiology

Boston Children's Hospital

Harvard Medical School



Boston Children's Hospital



HARVARD MEDICAL SCHOOL
TEACHING HOSPITAL



Declaration of Financial Interests or Relationships

Speaker Name: Simon K. Warfield, Ph.D.

I have no financial interests or relationships to disclose with regard to the subject matter of this presentation.



Motivation for Kidney DWI

- The kidney is a highly perfused organ with complex anatomy.
 - MRI can reveal morphology, vasculature and function.
- DWI is sensitive to kidney microstructure that changes with maturation, disease and injury:
 - CKD: tubular atrophy, tubular interstitial fibrosis
 - Renal ischemia
 - Acute pyelonephritis





Figure from Chapter 28 of "Principles of Biology" 2016 by Robert Bear and David Rintoul <http://cnx.org/content/m47445/1.4/>

- Goal: Improved imaging markers that better characterize microstructure
 - Caroli et al. Nephrol Dial Transplant (2018) doi:10.1093/ndt/gfy163




DWI for Kidney Imaging

- DWI contrast is manipulated by changing the diffusion encoding
 - Diffusion length scale, b-value
 - Orientation of diffusion gradient
- Microstructure from DWI
 - Comparison of signal changes reveals microstructural properties
- Motion can arise from
 - Respiration
 - Cardiac pulsation
 - Bulk motion
 - the subject rolls in the scanner.



DWI for Kidney Imaging: Motion

- Motion causes the signal from any one voxel in the images to arise from different places in the subject.
 - This reduces the fidelity of parameter estimates of models of the tissue microstructure.
- The DWI signal is measured at different times with varying gradient strengths and orientations.
 - We want to compare the DWI signals arising from the same anatomical location.
- Geometric distortion can also cause misalignment.




DWI for Kidney Imaging: Motion

- Approaches for motion compensation in renal DW-MRI:
- Navigator triggering, respiratory gating:
 - Acquire data at one instant of the respiratory cycle
 - Don't acquire data at other times in the respiratory cycle
 - Increases scan time (not very practical for routine clinical use)
- Breath holding:
 - Eliminate respiratory motion by not breathing during the scan.
 - Children and subjects who are ill may not be able to breath hold for long.
- We would prefer to continuously acquire DWI during free breathing:
 - Seek to minimize the total scan time
 - Increase patient comfort




DWI for Kidney Imaging

- Geometric distortion
 - Due to additional phase change during EPI readout due to magnetic susceptibility differences.



DWI for Kidney Imaging

- Geometric distortion
 - Due to additional phase change during EPI readout due to magnetic susceptibility differences.



- Image acquired with opposite phase encoding direction do not align.
 - Any change in orientation between the static magnetic field and the phase encoding direction alters the geometric distortion.
 - Magnitude of change in distortion depends on the local magnetic susceptibility differences.




Conventional Distortion Correction

- Measure the distortion field
- Apply correction for the distortion field
- Sequences:
 - PSF mapping
 - Gradient echo sequence to measure phase change
 - Spin echo with phase reversal.



Conventional Distortion Correction

- Measure the distortion field
- Apply correction for the distortion field
- Sequences:
 - PSF mapping
 - Gradient echo sequence to measure phase change
 - Spin echo with phase reversal
 - Acquire two EPI images
 - Opposite phase encoding directions
 - Align the distorted images to measure the magnitude of the change in distortion.
 - Estimate distortion field.



k-space

Coll-Font, J., Afacan, O., Hoge, S., Garg, H., Shashi, K., Marami, B., Gholipour, A., Chow, J., Warfield, S. and Kurugol, S., 2021. Retrospective Distortion and Motion Correction for Free-Breathing DW-MRI of the Kidneys Using Dual-Echo EPI and Slice-to-Volume Registration. *Journal of Magnetic Resonance Imaging*, 53(5), pp.1432-1443.




Problem for Distortion Correction

- Distortion fields are not static!
 - Breathing moves the tissues and changes the distortion field
 - Pockets of air in the bowel move over time
 - Time to acquire two DWI volumes with opposite phase encoding direction is too long.
 - By the time the second phase encode is acquired, the anatomy has moved to a new position.

Proposed solution:

Acquire two readouts at the same position with opposite phase encoding direction.

EPI readout L->R



Coll-Font, J., et. al. Journal of Magnetic Resonance Imaging, 53(5), pp.1432-1443.



Proposed Distortion Correction

- Modify sequence to obtain images with a dual-echo EPI readout¹
 - Reverse phase encoding (L->R and R->L) on second readout
 - Very short delay (~40ms) between readouts
 - **Assume no motion in between two echos => distortion field is static**
- Correct for distortion using L->R and R->L readouts
 - T2 decay between readouts may be non-negligible!
 - Adjust for contrast differences
 - Apply median filter to improve SNR
 - Compensate for amplitude difference
 - Run registration-based distortion compensation²

[1] Afacan et al. *J. NeuroIm* 2020. DOI: 10.1111/jon.12708

[2] Hedouin et al. "Block-Matching Distortion Correction of Echo-Planar Images With Opposite Phase Encoding Directions" *IEEE TMI* 2017. DOI: 10.1109/TMI.2016.2646920






Proposed Distortion Correction

L->R readout




T2 HASTE

R->L readout

Subject 1



Subject 2



Distorted

Aligned

Distorted



Boston Children's Hospital





HARVARD MEDICAL SCHOOL
TEACHING HOSPITAL



Motion Compensation

- Respiration creates inter-slice motion.
- Slice-to-Volume Registration:
 - Segment each kidney to target organ moving with six degrees of freedom.
 - Motion estimate with registration of each slice to 3D volume.
 - Rigid body motion of kidney with Kalman filter through time.



Kurugol, S. et al., 2017. Motion-robust spatially constrained parameter estimation in renal diffusion-weighted MRI by 3D motion tracking and correction of sequential slices. In *Molecular Imaging, Reconstruction and Analysis of Moving Body Organs, and Stroke Imaging and Treatment* (pp. 75-85). Springer, Cham.



Experiments


- 10 Subjects (ages 29-38, 3 females).
 - Healthy volunteers
 - Acquisition: Imaged with dual-echo DWI sequence
 - 17 directions each with 10 b-values each
 - Analysis pipeline
 - Geometric Distortion Correction of dual echo slices.
 - Slice to volume alignment for motion compensation.
- For IVIM and DTI models
 - Estimate parameters
 - Leave-one-direction-out cross-validation

IMAGING PARAMETERS:

- Dual-echo DWI sequence
- EPI acquisition
 - L->R and R->L phase encoding
- 10 b-values and 17 directions
- 3T Siemens Prisma
- TE1/TE2/TR 72ms/108ms/7000ms
- 18 coronal slices
- Res. 2.81x2.81x4 mm



Result: Distortion and Motion Correction



Correction for both the geometric distortion of each image and for the respiration induced motion of each image.

Coll-Font, J., et. al. Journal of Magnetic Resonance Imaging, 53(5), pp.1432-1443.




Boston Children's Hospital



HARVARD MEDICAL SCHOOL
TEACHING HOSPITAL



Result: Slice-to-Volume Registration



Coll-Font, J., et. al. Journal of Magnetic Resonance Imaging, 53(5), pp.1432-1443.



Boston Children's Hospital



HARVARD MEDICAL SCHOOL
TEACHING HOSPITAL



Results: Parameter Estimates


Model param CV (%) / method		No correction	DiMoCo
IVIM	Slow diff. (D) [%]	25.7+/-18.3	28.6+/-18.4
	Fast diff. (D*) [%]	63.3+/-21.6	49.6+/-20.5
	Perf. Frac, (f) [%]	38.5+/-14.5	31.4+/-13.7
	nRMSE	0.290+/-0.052	0.154+/-0.056
DTI	MD [%]	8.1+/-1.8	4.3+/-2.6
	FA [%]	14.4+/-1.1	12.4+/-2.4
	nRMSE	0.28+/-0.04	0.12+/-0.01

Correction for distortion and motion leads to:

- Reduction in coefficient of variation
- Improved normalized root mean square error.




Results: IVIM parameter maps



The corrected data leads to improved IVIM parameter estimates.



Results: DTI parameters



The corrected data leads to improved DTI parameter estimates.




Results: DTI parameters

FA

FA

FA

Orientation



No correction

Volume alignment.

Slice alignment improved orientation
visualization.



Boston Children's Hospital



HARVARD MEDICAL SCHOOL
TEACHING HOSPITAL



Conclusions and Future Work

- Conclusions
 - We can correct for distortion despite the presence of motion
 - Enabled by dual echo DWI with polarity reversal for each slice.
 - We can correct for respiration induced changes in position.
 - The distortion correction makes it easier to compensate for the respiratory motion.
 - These corrections enable DWI to be aligned, increasing the fidelity of parameter estimates for models of tissue microstructure.

- Future work
 - We plan to evaluate the utility of models of kidney microstructure.
 - Reduced influence of the confounds due to motion and distortion.



Acknowledgements

- Ideas, data, figures, algorithms, software created together with:
 - Sila Kurugol
 - Jaume Coll-font
 - Onur Afacan
 - Scott Hoge
 - Bahram Marami
 - Ali Gholipour
 - Jeanne Chow
- Supported in part by NIH awards:
 - R01 DK125561, R21 DK123569, R21 EB029627, R01 DK100404, R01 EB018988, R01 EB019483 and S10 OD025111.



Visit us at crl.med.harvard.edu



Boston Children's Hospital



**HARVARD MEDICAL SCHOOL
TEACHING HOSPITAL**

